



TemplinTM
preserving life

Blood Bank Refrigerator



Blood Bank Refrigerator

● Features

- Capacity of 60 bags to 875
- Microprocessor controlled temperature controller.
- Interior illumination for working
- 7 Days circular chart recorder
- CFC Free Compressor.
- Specially designed stainless steel
- Outer Glass

BLOOD BANK REFRIGERATORS meet the strict compliance requirements for storage of whole blood I pack cells needed in the blood banks, storage centers or hospitals. They are manufactured with highest precision in our modern manufacturing facility by the deft professionals using the high grade components and superior technology. It is available indifferent capacities and can also be customized.

Material of Construction

Inner stainless steel (304) 0.8mm thickness and outer mild steel, 1.0m m thickness powder coated.

Model	Capacity in Ltrs	Approx no. of Bags	No.of Baskets	Inner Dimensions H x W x D (In Cm)
TI-BR-10	100	60	2	60 x 43 x 40
TI-BR-11	200	145	4	80 x 50 x 50
TI-BR-12	300	240	4	100 x 58 x 50
TI-BR-13	400	300	4	105 x 60 x 60
TI-BR-14	500	350	5	120 x 60 x 65
TI-BR-15	600	425	5	125 x 68 x 71
TI-BR-19	1200	875	7 + 7	125 x 125 x 77



Email Us

info@tempinstruments.com
sales@tempinstruments.com

Address: 260 Salmon St,London NW9 8XY, United Kingdom