



TemplinTM
preserving life

Platelet Agitator



Platelet Agitator

● Features

- Available in 24 bags, 36 bags and 48 bags capacity.
- Movement of 70 ± 2 strokes/ minutes.
- Heavy duty sliders for trays and ball bearing.

The Platelet Agitator features the uniquely designed rack with high ventilation capability to allow oxygen flow to the platelets. Smooth sliding stainless steel trays are provided which can be removed for cleaning as well as for placing of the bags. Fixed speed agitation with heavy duty bearings and quiet operation also marked as special aspects in the design. The agitation stops automatically when the door is opened.

Technical Specification

| Model no. | TI- P - 125 | TI- P - 250 |
|----------------------|----------------------|----------------------|
| Capacity | 36 bags | 48 bags |
| Oscillation (stroke) | 70/Stroke Per Minute | 70/Stroke Per Minute |
| Power Supply | 230 VAC 50 Hz | 230 VAC 50 Hz |
| Alarm | Motion Failure | Motion Failure |



Email Us
info@tempinstruments.com
sales@tempinstruments.com

Address: 260 Salmon St, London NW9 8XY, United Kingdom