



TemplinTM
preserving life

TI-C-42



TI-C-42

● Features

- Precise Temperature Control and Safe Storage
- Visual & Audio Alarm system
- Dual Function – AC/DC
- Human-oriented
- Energy-saving & Environmentally Friendly
- Battery Safety

Specification

Model & Title	TI-C-42
Capacity in Litres	42
External Dim. W*D*H(mm)	364 x 680 x 515
Temperature. Range (°C)	10° ~ -20°
Cooling Mode	Direct cooling
Insulation Thickness	63
Climate Class	35° C
Weight	23.5 Kgs
Input Power	60W
Microprocessor Controller	Yes
Inside Light	Yes
Refrigeration	R134a
Alarm System - Door Open	Yes
Alarm System - Low Voltage	Yes
Power Supply (V/Hz) 12/24V	In-built
Power Supply (V/Hz) 220V	In-built



Email Us

info@tempinstruments.com

sales@tempinstruments.com

Address: 260 Salmon St, London NW9 8XY, United Kingdom