



**Templin**<sup>TM</sup>  
*preserving life*

**TI-C-60**



## TI-C-60

### ● Features

- Precise Temperature Control and Safe Storage
- Visual & Audio Alarm system
- Dual Function – AC/DC
- Human-oriented
- Energy-saving & Environmentally Friendly
- Battery Safety

## Specification

Model & Title	TI-C-60
Capacity in Litres	60
External Dim. W*D*H(mm)	740 x 495 x 520
Temperature. Range (°C)	10° ~ -20°
Cooling Mode	Direct cooling
Insulation Thickness	63
Climate Class	35° C
Weight	28 Kgs
Input Power	65 W
Microprocessor Controller	Yes
Temperature Display	DigitalDisplay
Inside Light	Yes
Refrigeration	R134a
Alarm System - Door Open	Yes
Alarm System - Low Voltage	Yes
Power Supply (V/Hz) 12/24V	In-built
Power Supply (V/Hz) 220V	In-built
External Lock (Latch)	Yes



Email Us

[info@tempinstruments.com](mailto:info@tempinstruments.com)

[sales@tempinstruments.com](mailto:sales@tempinstruments.com)

Address: 260 Salmon St, London NW9 8XY, United Kingdom