



**Templin**<sup>TM</sup>  
*preserving life*

**Thawing Bath (37°C)**



## Thawing Bath (37°C)

### Features

- Capacity of 8 bags & 15
- Microprocessor controlled temperature controller.
- CFC Free Compressor.
- Specially designed stainless steel

## Material of Construction

- Inner stainless steel (304)
- 0.8mm thickness and outer mild steel
- 1.0m m thickness powder coated.

## Technical Specification

Model no.	Capacity in liters	Approx no. of bags
TI-TB-8	10	8
TI-TB-15	20	15



Email Us  
[info@tempinstruments.com](mailto:info@tempinstruments.com)  
[sales@tempinstruments.com](mailto:sales@tempinstruments.com)

Address: 260 Salmon St, London NW9 8XY, United Kingdom