



TemplinTM
preserving life

Autoclave



Verticle high pressure autoclaves

- Fully Stainless Steel
 - Semi Automatic and fully automatic models
 - Double safety with safety valve & Automatic digital controller
 - Automatic low waterlevel cut-off device to protect heaters
 - Single, Double & Triple walled (jacked)
- autoclaves are extremely compact, sturdy and decent in look

Features

- Top end model for sterilization of all critical application with exact reproducibility
- Semi/ fully automatic operation
- Unique double / tripped wall design
- Heavy duty jacket
- Automatic purging of steam
- Micro processor based precise controlling
- Set pressure & selected temperature are attained
- Timer starts countdown after attaining set temperature
- Audio alarm buzzer gives signal of completion of sterilization cycle after set time is over
- Recording of one cycle - optional

Specification

Model	TI-EV-15	TI-EV-24	TI-EV-35	TI-EV-41
Size - Internal	12" Dia x 20" Depth	16" Dia x 24" Depth	20" Dia x 28" Depth	22" Dia x 30" Depth
Construction of Chamber	S.S 304	S.S 304	S.S 304	S.S 304
Capacity	35 Litres	75 Litres	150 Litres	175 Litres
Power consumption	3000 Watts	4000 Watts	6000 Watts	9000 Watts
operating pressure	1.2 kg	1.2 kg	1.2 kg	1.2 kg
operating temperature	121 C	121 C	121 C	121 C
voltage	230 volt	230 volt	440 volt	440 volt
control	Digital	Digital	Digital	Digital

Tempin[™]
preserving life

Email Us

info@tempinstruments.com
sales@tempinstruments.com

Address: 260 Salmon St, London NW9 8XY, United Kingdom